

Rcn2 Cas9-KO Strategy

Designer:

Daohua Xu

Reviewer:

Huimin Su

Design Date:

2020-1-13

Project Overview

Project Name

Rcn2

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Rcn2* gene. The schematic diagram is as follows:



- The *Rcn2* gene has 4 transcripts. According to the structure of *Rcn2* gene, exon3-exon4 of *Rcn2-201* (ENSMUST00000114276.2) transcript is recommended as the knockout region. The region contains 311bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Rcn2* gene. The brief process is as follows: CRISPR/Cas9 system v

- The *Rcn2* gene is located on the Chr9. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Rcn2 reticulocalbin 2 [Mus musculus (house mouse)]

Gene ID: 26611, updated on 31-Jan-2019

Summary




Official Symbol	Rcn2 provided by MGI
Official Full Name	reticulocalbin 2 provided by MGI
Primary source	MGI:MGI:1349765
See related	Ensembl:ENSMUSG00000032320
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	AA408742, Tcbp49
Expression	Ubiquitous expression in CNS E11.5 (RPKM 45.9), limb E14.5 (RPKM 44.7) and 27 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Rcn2-201	ENSMUST00000114276.2	2024	320aa	Protein coding	CCDS23206	Q8BP92	TSL:1 GENCODE basic APPRIS P1
Rcn2-203	ENSMUST00000147842.1	3080	71aa	Nonsense mediated decay	-	D6RHL9	TSL:1
Rcn2-204	ENSMUST00000151585.1	1987	No protein	Retained intron	-	-	TSL:1
Rcn2-202	ENSMUST00000144869.1	655	No protein	Retained intron	-	-	TSL:3


The strategy is based on the design of *Rcn2-201* transcript,The transcription is shown below



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

