

# ***Abhd11*** Cas9-KO Strategy

Designer: Xueting Zhang

Reviewer: Yanhua Shen

Date: 2020-03-11

# Project Overview

**Project Name**

*Abhd11*

**Project type**


**Cas9-KO**

**Strain background**

**C57BL/6JGpt**

# Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Abhd11* gene. The schematic diagram is as follows:



- The *Abhd11* gene has 7 transcripts. According to the structure of *Abhd11* gene, exon3-exon6 of *Abhd11-201* (ENSMUST00000046999.11) transcript is recommended as the knockout region. The region contains most of the coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Abhd11* gene. The brief process is as follows: CRISPR/Cas9 system

- *Abhd11os* gene will be destroyed.
- The knockout region is near to the N-terminal of *Gm43372* gene, this strategy may influence the regulatory function of the N-terminal of *Gm43372* gene.
- The *Abhd11* gene is located on the Chr5. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This Strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

# Gene information (NCBI)

## Abhd11 abhydrolase domain containing 11 [ *Mus musculus* (house mouse) ]

Gene ID: 68758, updated on 27-Jan-2020

### Summary

**Official Symbol** Abhd11 provided by [MGI](#)  
**Official Full Name** abhydrolase domain containing 11 provided by [MGI](#)  
**Primary source** [MGI:MGI:1916008](#)  
**See related** [Ensembl:ENSMUSG00000040532](#)  
**Gene type** protein coding  
**RefSeq status** VALIDATED  
**Organism** [Mus musculus](#)  
**Lineage** Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus  
**Also known as** Wbscr21; 1110054D16Rik; A630008N09Rik  
**Expression** Ubiquitous expression in large intestine adult (RPKM 44.7), colon adult (RPKM 41.1) and 28 other tissues [See more](#)  
**Orthologs** [human](#) [all](#)

### Genomic context

**Location:** 5; 5 G2

See Abhd11 in [Genome Data Viewer](#)

**Exon count:** 6


| Annotation release  | Status            | Assembly                                       | Chr | Location                           |
|---------------------|-------------------|------------------------------------------------|-----|------------------------------------|
| <a href="#">108</a> | current           | GRCm38.p6 ( <a href="#">GCF_000001635.26</a> ) | 5   | NC_000071.6 (135009152..135012175) |
| Build 37.2          | previous assembly | MGSCv37 ( <a href="#">GCF_000001635.18</a> )   | 5   | NC_000071.5 (135485022..135488045) |

# Transcript information (Ensembl)


The gene has 7 transcripts,all transcripts are shown below:

| Name       | Transcript ID                         | bp   | Protein               | Biotype                 | CCDS                      | UniProt                    | Flags                         |
|------------|---------------------------------------|------|-----------------------|-------------------------|---------------------------|----------------------------|-------------------------------|
| Abhd11-201 | <a href="#">ENSMUST00000046999.11</a> | 1493 | <a href="#">307aa</a> | Protein coding          | <a href="#">CCDS19730</a> | <a href="#">Q8K4F5</a>     | TSL:1 GENCODE basic APPRIS P1 |
| Abhd11-206 | <a href="#">ENSMUST00000154469.7</a>  | 1473 | <a href="#">227aa</a> | Protein coding          | <a href="#">CCDS51659</a> | <a href="#">D3YYK0</a>     | TSL:1 GENCODE basic           |
| Abhd11-202 | <a href="#">ENSMUST00000111216.7</a>  | 972  | <a href="#">227aa</a> | Protein coding          | <a href="#">CCDS51659</a> | <a href="#">D3YYK0</a>     | TSL:3 GENCODE basic           |
| Abhd11-207 | <a href="#">ENSMUST00000201890.1</a>  | 476  | <a href="#">59aa</a>  | Protein coding          | -                         | <a href="#">A0A0J9YUT1</a> | CDS 3' incomplete TSL:5       |
| Abhd11-205 | <a href="#">ENSMUST00000149778.1</a>  | 1475 | <a href="#">78aa</a>  | Nonsense mediated decay | -                         | <a href="#">A0A0J9YV94</a> | TSL:2                         |
| Abhd11-203 | <a href="#">ENSMUST00000123227.2</a>  | 1466 | <a href="#">158aa</a> | Nonsense mediated decay | -                         | <a href="#">A0A0J9YU40</a> | TSL:5                         |
| Abhd11-204 | <a href="#">ENSMUST00000148831.2</a>  | 731  | <a href="#">67aa</a>  | Nonsense mediated decay | -                         | <a href="#">A0A0J9YVC2</a> | TSL:3                         |


The strategy is based on the design of *Abhd11-201* transcript,The transcription is shown below




# Genomic location distribution



# Protein domain



# Mouse phenotype description(MGI )



*Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).*

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

