

Calml3 Cas9-KO Strategy

Designer: Daohua Xu

Reviewer: Huimin Su

Design Date: 2020-5-26

Project Overview

Project Name

Calml3

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Calml3* gene. The schematic diagram is as follows:



- The *Calml3* gene has 2 transcripts. According to the structure of *Calml3* gene, exon1 of *Calml3-201* (ENSMUST00000077698.4) transcript is recommended as the knockout region. The region contains all of the coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Calml3* gene. The brief process is as follows: CRISPR/Cas9 system

- The *Calml3* gene is located on the Chr13. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Calml3 calmodulin-like 3 [Mus musculus (house mouse)]

Gene ID: 70405, updated on 13-Mar-2020

Summary




Official Symbol	Calml3 provided by MGI
Official Full Name	calmodulin-like 3 provided by MGI
Primary source	MGI:MGI:1917655
See related	Ensembl:ENSMUSG00000063130
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	2310068O22Rik, A1326174
Orthologs	human all

Transcript information (Ensembl)


The gene has 2 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Calml3-201	ENSMUST00000077698.4	1421	149aa	Protein coding	CCDS26214	Q9D6P8	TSL:NA GENCODE basic APPRIS P1
Calml3-202	ENSMUST00000222976.1	509	115aa	Protein coding	-	A0A1Y7VK66	CDS 5' incomplete TSL:5


The strategy is based on the design of *Calml3-201* transcript,the transcription is shown below:



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

