

Unc45b Cas9-KO Strategy

Designer: JiaYu

Reviewer: Xiaojing Li

Design Date: 2020-7-29

Project Overview



Project Name

Unc45b

Project type

Cas9-KO


Strain background

C57BL/6JGpt

Knockout strategy



This model will use CRISPR/Cas9 technology to edit the *Unc45b* gene. The schematic diagram is as follows:



Technical routes



- ➤ The *Unc45b* gene has 5 transcripts. According to the structure of *Unc45b* gene, exon4 of *Unc45b*201(ENSMUST00000018989.13) transcript is recommended as the knockout region. The region contains 176bp coding sequence. Knock out the region will result in disruption of protein function.
- ➤ In this project we use CRISPR/Cas9 technology to modify *Unc45b* gene. The brief process is as follows: CRISPR/Cas9 system were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

Notice



- > According to the existing MGI data, mice homozygous for a transgenic gene disruption exhibit embryonic lethality at E9 without placental abnormalities.
- The *Unc45b* gene is located on the Chr11. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- > This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)



Unc45b unc-45 myosin chaperone B [Mus musculus (house mouse)]

Gene ID: 217012, updated on 13-Mar-2020

Summary

☆ ?

Official Symbol Unc45b provided by MGI

Official Full Name unc-45 myosin chaperone B provided by MGI

Primary source MGI:MGI:2443377

See related Ensembl:ENSMUSG00000018845

Gene type protein coding
RefSeq status VALIDATED
Organism Mus musculus

Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia;

Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus

Also known as AA445617, Cmya4, D230041A13Rik, Unc45

Expression Biased expression in heart adult (RPKM 41.6), mammary gland adult (RPKM 6.5) and 5 other tissuesSee more

Orthologs <u>human all</u>


Transcript information (Ensembl)



The gene has 5 transcripts, all transcripts are shown below:


Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Unc45b-201	ENSMUST00000018989.13	3549	929aa	Protein coding	CCDS36251	Q8CGY6	TSL:1 GENCODE basic APPRIS P2
Unc45b-202	ENSMUST00000108160.7	4152	<u>931aa</u>	Protein coding	÷	Q8CGY6	TSL:5 GENCODE basic APPRIS ALT1
Unc45b-205	ENSMUST00000164945.2	3555	<u>931aa</u>	Protein coding	20	Q8CGY6	TSL:5 GENCODE basic APPRIS ALT1
Unc45b-204	ENSMUST00000142336.1	498	No protein	Processed transcript	51	-	TSL:5
Unc45b-203	ENSMUST00000125499.1	446	No protein	Processed transcript	21	-	TSL:5

The strategy is based on the design of *Unc45b-201* transcript, the transcription is shown below:




Genomic location distribution






Protein domain





Mouse phenotype description(MGI)





Phenotypes affected by the gene are marked in blue.Data quoted from MGI database(http://www.informatics.jax.org/).

According to the existing MGI data, mice homozygous for a transgenic gene disruption exhibit embryonic lethality at E9 without placental abnormalities.



If you have any questions, you are welcome to inquire. Tel: 400-9660890





