

*Efh***d1** Cas9-KO Strategy

Designer: JiaYu

Reviewer: Xiaojing Li

Design Date: 2020-8-27

Project Overview

Project Name

Efhd1

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Efhdl* gene. The schematic diagram is as follows:



- The *Efhdl* gene has 3 transcripts. According to the structure of *Efhdl* gene, exon2 of *Efhdl*-201(ENSMUST00000027472.6) transcript is recommended as the knockout region. The region contains 148bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Efhdl* gene. The brief process is as follows: CRISPR/Cas9 system were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

- The *Efhdl* gene is located on the Chr1. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Efhd1 EF hand domain containing 1 [Mus musculus (house mouse)]

Gene ID: 98363, updated on 13-Mar-2020

Summary




Official Symbol	Efhd1 provided by MGI
Official Full Name	EF hand domain containing 1 provided by MGI
Primary source	MGI:MGI:1921607
See related	Ensembl:ENSMUSG00000026255
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	4931430I01Rik, AI452351, PP3051, mitocalcin
Expression	Biased expression in adrenal adult (RPKM 265.2), testis adult (RPKM 183.7) and 4 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 3 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Efhd1-201	ENSMUST00000027472.6	1854	240aa	Protein coding	CCDS15132	Q9D4J1	TSL:1 GENCODE basic APPRIS P1
Efhd1-202	ENSMUST00000118687.7	1726	240aa	Nonsense mediated decay	-	Q9D4J1	TSL:1
Efhd1-203	ENSMUST00000150831.1	353	No protein	Processed transcript	-	-	TSL:2


The strategy is based on the design of *Efhd1-201* transcript,the transcription is shown below:



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

