

# ***Samd7 Cas9-KO Strategy***

**Designer: Zihe Cui**

**Reviewer: Yanhua Shen**

**Design Date: 2020-9-30**

# Project Overview

**Project Name**

***Samd7***

**Project type**


**Cas9-KO**

**Strain background**

**C57BL/6JGpt**

# Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Samd7* gene. The schematic diagram is as follows:



- The *Samd7* gene has 4 transcripts. According to the structure of *Samd7* gene, exon3-exon6 of *Samd7*-201(ENSMUST00000108262.9) transcript is recommended as the knockout region. The region contains 934bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Samd7* gene. The brief process is as follows: CRISPR/Cas9 system were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

- According to the existing MGI data, homozygous knockout reduces the sensitivity of the retinal rods to low-to-moderate flash luminescence.
- The *Samd7* gene is located on the Chr3. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.



# Gene information (NCBI)

## Samd7 sterile alpha motif domain containing 7 [ *Mus musculus* (house mouse) ]

Gene ID: 75953, updated on 25-Sep-2020

### Summary


Official Symbol	Samd7 provided by <a href="#">MGI</a>
Official Full Name	sterile alpha motif domain containing 7 provided by <a href="#">MGI</a>
Primary source	<a href="#">MGI:MGI:1923203</a>
See related	<a href="#">Ensembl:ENSMUSG00000051860</a>
Gene type	protein coding
RefSeq status	VALIDATED
Organism	<a href="#">Mus musculus</a>
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	4930597A01Rik
Expression	Biased expression in testis adult (RPKM 1.3) and CNS E18 (RPKM 0.1) <a href="#">See more</a>
Orthologs	<a href="#">human</a> <a href="#">all</a>

# Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Samd7-201	<a href="#">ENSMUST00000108262.9</a>	3144	<a href="#">445aa</a>	Protein coding	<a href="#">CCDS17285</a>	<a href="#">Q8C8Y5</a>	TSL:5 GENCODE basic APPRIS P1
Samd7-203	<a href="#">ENSMUST00000174395.7</a>	3058	<a href="#">39aa</a>	Nonsense mediated decay	-	<a href="#">Q8C8Y5</a>	TSL:2
Samd7-202	<a href="#">ENSMUST00000172593.1</a>	1511	<a href="#">102aa</a>	Nonsense mediated decay	-	<a href="#">G3UWL1</a>	TSL:5
Samd7-204	<a href="#">ENSMUST00000192872.1</a>	1241	No protein	Retained intron	-	-	TSL:NA

The strategy is based on the design of *Samd7-201* transcript,the transcription is shown below:




# Genomic location distribution






# Protein domain



# Mouse phenotype description(MGI )



*Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).*

According to the existing MGI data, homozygous knockout reduces the sensitivity of the retinal rods to low-to-moderate flash luminescence.

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

