

Marchf6 Cas9-KO Strategy

Designer: Xueting Zhang

Reviewer: Daohua Xu

Design Date: 2020-9-16

Project Overview

Project Name

Marchf6

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Marchf6* gene. The schematic diagram is as follows:



- The *Marchf6* gene has 4 transcripts. According to the structure of *Marchf6* gene, exon2 of *March6-201*(ENSMUST00000090227.5) transcript is recommended as the knockout region. The region contains 97bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Marchf6* gene. The brief process is as follows: CRISPR/Cas9 system were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

- The *Marchf6* gene is located on the Chr15. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- Transcript *Marchf6*-202&203 may not be affected.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Gene information (NCBI)

Marchf6 membrane associated ring-CH-type finger 6 [Mus musculus (house mouse)]

Gene ID: 223455, updated on 13-Mar-2020

Summary




Official Symbol	Marchf6 provided by MGI
Official Full Name	membrane associated ring-CH-type finger 6 provided by MGI
Primary source	MGI:MGI:2442773
See related	Ensembl:ENSMUSG00000039100
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	3830408G03, F830029L24Rik, March6, mKIAA0597
Expression	Ubiquitous expression in CNS E18 (RPKM 28.9), whole brain E14.5 (RPKM 26.6) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
March6-201	ENSMUST00000090227.5	6260	909aa	Protein coding	CCDS37054	Q6ZQ89	TSL:1 GENCODE basic APPRIS P1
March6-204	ENSMUST00000228122.1	662	No protein	Processed transcript	-	-	
March6-202	ENSMUST00000226207.1	2669	No protein	Retained intron	-	-	
March6-203	ENSMUST00000226868.1	736	No protein	Retained intron	-	-	


The strategy is based on the design of *March6-201* transcript,the transcription is shown below:



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

