

Cpn2 Cas9-KO Strategy

Designer: Zihe Cui

Reviewer: Xueting Zhang

Design Date: 2020-10-10

Project Overview

Project Name

Cpn2

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Cpn2* gene. The schematic diagram is as follows:



- The *Cpn2* gene has 1 transcript. According to the structure of *Cpn2* gene, exon2 of *Cpn2-201*(ENSMUST00000064856.8) transcript is recommended as the knockout region. The region contains all of the coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Cpn2* gene. The brief process is as follows: CRISPR/Cas9 system were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

- The KO region is close to *Gm49645* gene. Knockout the region may affect the function of *Gm49645* gene.
- The *Cpn2* gene is located on the Chr16. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.



Gene information (NCBI)

Cpn2 carboxypeptidase N, polypeptide 2 [*Mus musculus* (house mouse)]

Gene ID: 71756, updated on 25-Sep-2020

Summary




Official Symbol	Cpn2 provided by MGI
Official Full Name	carboxypeptidase N, polypeptide 2 provided by MGI
Primary source	MGI:MGI:1919006
See related	Ensembl:ENSMUSG00000023176
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	1300018K11Rik
Expression	Biased expression in liver adult (RPKM 111.2) and liver E18 (RPKM 15.0) See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 1 transcript, and the transcript is shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Cpn2-201	ENSMUST00000064856.8	4550	547aa	Protein coding	CCDS49819	Q9DBB9	TSL:1 GENCODE basic APPRIS P1


The strategy is based on the design of *Cpn2-201* transcript, the transcription is shown below:



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.

Tel: 400-9660890

