

Chmp7 Cas9-KO Strategy

Designer: Xueting Zhang

Reviewer: Daohua Xu

Design Date: 2020-9-29

Project Overview

Project Name

Chmp7

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Chmp7* gene. The schematic diagram is as follows:



- The *Chmp7* gene has 4 transcripts. According to the structure of *Chmp7* gene, exon2 of *Chmp7*-201(ENSMUST00000036381.10) transcript is recommended as the knockout region. The region contains 172bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Chmp7* gene. The brief process is as follows: CRISPR/Cas9 system were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

- The *Chmp7* gene is located on the Chr14. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Chmp7 charged multivesicular body protein 7 [Mus musculus (house mouse)]

Gene ID: 105513, updated on 13-Mar-2020

Summary




Official Symbol	Chmp7 provided by MGI
Official Full Name	charged multivesicular body protein 7 provided by MGI
Primary source	MGI:MGI:1913922
See related	Ensembl:ENSMUSG00000034190
Gene type	protein coding
RefSeq status	VALIDATED
Organism	Mus musculus
Lineage	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Also known as	4930596K11Rik, 6330407G04Rik, AI450338, AW550775
Expression	Ubiquitous expression in CNS E18 (RPKM 31.1), whole brain E14.5 (RPKM 26.1) and 28 other tissues See more
Orthologs	human all

Transcript information (Ensembl)


The gene has 4 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Chmp7-201	ENSMUST00000036381.9	2649	451aa	Protein coding	CCDS27242	Q8R1T1	TSL:1 GENCODE basic APPRIS P1
Chmp7-203	ENSMUST00000224964.1	3348	No protein	Retained intron	-	-	
Chmp7-204	ENSMUST00000225036.1	3115	No protein	Retained intron	-	-	
Chmp7-202	ENSMUST00000224601.1	2447	No protein	Retained intron	-	-	


The strategy is based on the design of *Chmp7-201* transcript,the transcription is shown below:




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue.Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.

Tel: 400-9660890

