

Ccl27a Cas9-KO Strategy

Designer: Huan Wang

Reviewer: Yumeng Wang

Design Date: 2021-9-26

Project Overview

Project Name

Ccl27a

Project type


Cas9-KO

Strain background

C57BL/6JGpt

Knockout strategy

This model will use CRISPR/Cas9 technology to edit the *Ccl27a* gene. The schematic diagram is as follows:



- The *Ccl27a* gene has 19 transcripts. According to the structure of *Ccl27a* gene, exon1-exon4 of *Ccl27a*-206(ENSMUST00000108037.9) transcript is recommended as the knockout region. The region contains 205bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR/Cas9 technology to modify *Ccl27a* gene. The brief process is as follows: CRISPR/Cas9 system were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.

- The *Ccl27a* gene is located on the Chr4. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of the gene knockout on gene transcription, RNA splicing and protein translation cannot be predicted at the existing technology level.

Ccl27a chemokine (C-C motif) ligand 27A [Mus musculus (house mouse)]

Gene ID: 20301, updated on 13-Mar-2020

Summary



Official Symbol Ccl27a provided by [MGI](#)

Official Full Name chemokine (C-C motif) ligand 27A provided by [MGI](#)

Primary source [MGI:MGI:1343459](#)

See related [Ensembl:ENSMUSG00000073888](#)

Gene type protein coding

RefSeq status VALIDATED

Organism [Mus musculus](#)

Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus

Also known as ALP, AW558992, CTACK, CTAK, Ccl27, ESkin, ILC, PESKY, Scya27, Scya27a

Expression Broad expression in testis adult (RPKM 85.9), cortex adult (RPKM 62.4) and 21 other tissues [See more](#)


Orthologs [human](#) [all](#)

Transcript information（Ensembl）


The gene has 19 transcripts,all transcripts are shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Ccl27a-205	ENSMUST00000108036.7	841	164aa	Protein coding	CCDS38721	A2AMS7	TSL:2 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS ALT2
Ccl27a-204	ENSMUST00000108035.7	731	127aa	Protein coding	CCDS84714	A2AMS6	TSL:1 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS ALT2
Ccl27a-219	ENSMUST00000173865.8	657	87aa	Protein coding	CCDS51144	A2AMS8	TSL:1 GENCODE basic
Ccl27a-206	ENSMUST00000108037.8	590	164aa	Protein coding	CCDS38721	A2AMS7	TSL:2 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS ALT2
Ccl27a-201	ENSMUST0000074387.3	405	120aa	Protein coding	CCDS18072	A0A0R4J0X6	TSL:1 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS P3
Ccl27a-202	ENSMUST00000108032.2	350	87aa	Protein coding	CCDS51144	A2AMS8	TSL:1 GENCODE basic
Ccl27a-217	ENSMUST00000155240.1	343	87aa	Protein coding	CCDS51144	A2AMS8	TSL:2 GENCODE basic
Ccl27a-203	ENSMUST00000108033.7	428	119aa	Protein coding	-	A2AMS5	TSL:1 GENCODE basic APPRIS is a system to annotate alternatively spliced transcripts based on a range of computational methods to identify the most functionally important transcript(s) of a gene. APPRIS ALT2
Ccl27a-215	ENSMUST00000151142.7	822	37aa	Nonsense mediated decay	-	A0A2I3BRY3	TSL:3
Ccl27a-208	ENSMUST00000131349.7	489	40aa	Nonsense mediated decay	-	A0A2I3BQE5	TSL:3
Ccl27a-218	ENSMUST00000155322.7	460	40aa	Nonsense mediated decay	-	A0A2I3BQE5	TSL:3
Ccl27a-214	ENSMUST00000144090.7	996	No protein	Retained intron	-	-	TSL:2
Ccl27a-210	ENSMUST00000133572.7	798	No protein	Retained intron	-	-	TSL:2
Ccl27a-216	ENSMUST00000151909.7	614	No protein	Retained intron	-	-	TSL:1
Ccl27a-211	ENSMUST00000136768.7	556	No protein	Retained intron	-	-	TSL:1
Ccl27a-212	ENSMUST00000138337.7	552	No protein	Retained intron	-	-	TSL:2
Ccl27a-213	ENSMUST00000143676.1	487	No protein	Retained intron	-	-	TSL:2
Ccl27a-209	ENSMUST00000133556.7	453	No protein	Retained intron	-	-	TSL:2
Ccl27a-207	ENSMUST00000126759.1	361	No protein	Retained intron	-	-	TSL:3


The strategy is based on the design of *Ccl27a-206* transcript,the transcription is shown below:



Genomic location distribution



Protein domain



If you have any questions, you are welcome to inquire.
Tel: 400-9660890

