

Pdzd8 Cas9-CKO Strategy

Designer: Xiaojing Li

Reviewer: Jia Yu

Design Date: 2022-5-26

Project Overview

Project Name

Pdzd8

Project type


Cas9-CKO

Strain background

C57BL/6JGpt

Conditional Knockout strategy

This model will use CRISPR-Cas9 technology to edit the *Pdzd8* gene. The schematic diagram is as follows:



- The *Pdzd8* gene has 1 transcript. According to the structure of *Pdzd8* gene, exon3 of *Pdzd8-201*(ENSMUST00000099274.4) transcript is recommended as the knockout region. The region contains 103bp coding sequence. Knock out the region will result in disruption of protein function.
- In this project we use CRISPR-Cas9 technology to modify *Pdzd8* gene. The brief process is as follows: CRISPR-Cas9 system and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice. Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- The flox mice will be knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

- The *Pdzd8* gene is located on the Chr19. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)

Pdzd8 PDZ domain containing 8 [Mus musculus (house mouse)]

Gene ID: 107368, updated on 24-Apr-2022

Summary



Official Symbol Pdzd8 provided by [MGI](#)

Official Full Name PDZ domain containing 8 provided by [MGI](#)

Primary source [MGI:MGI:2677270](#)

See related [Ensembl:ENSMUSG00000074746](#)

Gene type protein coding

RefSeq status VALIDATED

Organism [Mus musculus](#)

Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus

Also known as A630041P07Rik, AV329468, AW228944, Pdzk8

Expression Ubiquitous expression in testis adult (RPKM 15.4), ovary adult (RPKM 6.8) and 28 other tissues [See more](#)


Orthologs [human](#) [all](#)

Transcript information (Ensembl)


The gene has 1 transcript,and the transcript is shown below:

Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Pdzd8-201	ENSMUST00000099274.4	5982	1147aa	Protein coding	CCDS29936		TSL:1 , GENCODE basic , APPRIS P1 ,

The strategy is based on the design of *Pdzd8-201* transcript,the transcription is shown below:




Genomic location distribution



Protein domain



Mouse phenotype description(MGI)



Phenotypes affected by the gene are marked in blue. Data quoted from MGI database(<http://www.informatics.jax.org/>).

If you have any questions, you are welcome to inquire.
Tel: 400-9660890

