

Genotyping Report

Strain ID	T008299	Strain Type	CKO(Cas9)	Genetic Background	C57BL/6J
Designer	Tian'e Wang	Gene Name	<i>Wtap</i>		

1. Strategy of Genotyping



 LoxP  Primer ① ② PCR Reaction Number  Sequencing Region


Wild type: ①PCR reaction obtains a single WT band; ②PCR reaction obtains a single WT band.
 Heterozygote: ①PCR reaction obtains a WT band and a Targeted band; ②PCR reaction obtains a WT band and a Targeted band.
 Homozygote: ①PCR reaction obtains a single Targeted band; ②PCR reaction obtains a single Targeted band.
 Note: The sizes of WT and Targeted band are shown below.

2. Primer Information

PCR No.	Primer No.	Primer Name	Sequence	Band Size
①(5'arm)	F1	02278-Wtap-5wt-tF2	AGTGCTGGGATAACAGGCAC	WT:302bp Targeted:407bp
	R1	02278-Wtap-5wt-tR2	TGTTTGAAACAATTTCAACTTAGG	
②(3'arm)	F2	02278-Wtap-3wt-tF2	AACAGGGTATATTCAAGTTGGCC	WT:355bp Targeted:456bp
	R2	02278-Wtap-3wt-tR2	ATAGGTTGGGCATAGGTTACCTAC	

3. Gel Image & Conclusion

--



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ （-0.5℃/cycle）	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2019-12-23