


## Genotyping Report

|           |           |             |              |                    |             |
|-----------|-----------|-------------|--------------|--------------------|-------------|
| Strain ID | T009476   | Strain Type | CKO(Cas9)    | Genetic Background | C57BL/6JGpt |
| Designer  | Ya'nan Xu | Gene Name   | <i>Csf1r</i> |                    |             |

### 1. Strategy of Genotyping



Wild type: ①PCR reaction obtains a single WT band; ②PCR reaction obtains a single WT band.

Heterozygote: ①PCR reaction obtains a WT band and a Targeted band; ②PCR reaction obtains a WT band and a Targeted band.



Homozygote: ①PCR reaction obtains a single Targeted band; ②PCR reaction obtains a single Targeted band.

Note: The sizes of WT and Targeted band are shown below.

### 2. Primer Information

| PCR No.  | Primer No. | Primer Name           | Sequence                  | Band Size                  |
|----------|------------|-----------------------|---------------------------|----------------------------|
| ①(5'arm) | F1         | JS03336-Csf1r-5wt-tF1 | AGCCTCCTTCCCTAAGACCTGACA  | WT:352bp<br>Targeted:454bp |
|          | R1         | JS03336-Csf1r-5wt-tR1 | TACCACTGGGCATCGGGTCTGT    |                            |
| ②(3'arm) | F2         | JS03336-Csf1r-3wt-tF1 | CCCTTCTGACAGGAAGAGAAGTTCC | WT:281bp<br>Targeted:384bp |
|          | R2         | JS03336-Csf1r-3wt-tR1 | CGAGGGACTGAACTCTGAGTCTCA  |                            |

### 3. Gel Image & Conclusion



Note: P: Positive control; WT: Wildtype control; B: Blank control (ddH<sub>2</sub>O); M: DNA Ladder

#### 4. PCR Condition

| PCR Reaction Component |                                                 |      |             |
|------------------------|-------------------------------------------------|------|-------------|
| Seg.                   | reaction component                              |      | Volume (μl) |
| 1                      | 2 × Taq Master Mix , Dye Plus, (Vazyme P112-03) |      | 12.5        |
| 2                      | ddH2O                                           |      | 9.5         |
| 3                      | Primer A(10pmol/μl)                             |      | 1           |
| 4                      | Primer B(10pmol/μl)                             |      | 1           |
| 5                      | Template(≈100ng/μl)                             |      | 1           |
| PCR program            |                                                 |      |             |
| Seg.                   | Temp.                                           | Time | Cycle       |
| 1                      | 95℃                                             | 5min | 20×         |
| 2                      | 98℃                                             | 30s  |             |
| 3                      | 65℃ （-0.5℃/cycle）                               | 30s  |             |
| 4                      | 72℃                                             | 45s  |             |
| 5                      | 98℃                                             | 30s  | 20×         |
| 6                      | 55℃                                             | 30s  |             |
| 7                      | 72℃                                             | 45s  |             |

|   |     |      |  |
|---|-----|------|--|
| 8 | 72℃ | 5min |  |
| 9 | 10℃ | hold |  |

Reviewer:Ting Sun

Date:2020-5-12