

## Genotyping Report

|           |           |             |           |                    |             |
|-----------|-----------|-------------|-----------|--------------------|-------------|
| Strain ID | T009951   | Strain Type | CKO(Cas9) | Genetic Background | C57BL/6JGpt |
| Designer  | Zifan Lin | Gene Name   | Enpp6     |                    |             |

### 1. Strategy of Genotyping



Wild type: ①PCR reaction obtains a single WT band; ②PCR reaction obtains a single WT band.

Heterozygote: ①PCR reaction obtains a WT band and a Targeted band; ②PCR reaction obtains a WT band and a Targeted band.


Homozygote: ①PCR reaction obtains a single Targeted band; ②PCR reaction obtains a single Targeted band.

Note: The sizes of WT and Targeted band are shown below.

### 2. Primer Information

| PCR No.  | Primer No. | Primer Name           | Sequence                 | Band Size                  |
|----------|------------|-----------------------|--------------------------|----------------------------|
| ①(5'arm) | F1         | JS03615-Enpp6-5wt-tF1 | TTCCCTGAACTGTAGGCAAGAGC  | WT:386bp<br>Targeted:491bp |
|          | R1         | JS03615-Enpp6-5wt-tR1 | AGCCAGTGTCTCCTTCTCAGAGC  |                            |
| ②(3'arm) | F2         | JS03615-Enpp6-3wt-tF1 | CTTGTTTGAAATGTGCAAAGCAC  | WT:343bp<br>Targeted:444bp |
|          | R2         | JS03615-Enpp6-3wt-tR1 | ATCTGTGTCTTCTGCTTTGACTCC |                            |

### 3. Gel Image & Conclusion



Note: P: Positive control; WT: Wildtype control; B: Blank control (ddH<sub>2</sub>O); M: DNA Ladder

#### 4. PCR Condition

| PCR Reaction Component |                                                 |      |             |
|------------------------|-------------------------------------------------|------|-------------|
| Seg.                   | reaction component                              |      | Volume (μl) |
| 1                      | 2 × Taq Master Mix , Dye Plus, (Vazyme P112-03) |      | 12.5        |
| 2                      | ddH2O                                           |      | 9.5         |
| 3                      | Primer A(10pmol/μl)                             |      | 1           |
| 4                      | Primer B(10pmol/μl)                             |      | 1           |
| 5                      | Template(≈100ng/μl)                             |      | 1           |
| PCR program            |                                                 |      |             |
| Seg.                   | Temp.                                           | Time | Cycle       |
| 1                      | 95℃                                             | 5min | 20×         |
| 2                      | 98℃                                             | 30s  |             |
| 3                      | 65℃ （-0.5℃/cycle）                               | 30s  |             |
| 4                      | 72℃                                             | 45s  |             |
| 5                      | 98℃                                             | 30s  | 20×         |
| 6                      | 55℃                                             | 30s  |             |
| 7                      | 72℃                                             | 45s  |             |
| 8                      | 72℃                                             | 5min |             |
| 9                      | 10℃                                             | hold |             |

Reviewer:Ting Sun

Date:2020-5-12