

Genotyping Report

Strain ID	T018255	Strain Type	CKO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Ya'nan Xu	Gene Name	Homer1		

1. Strategy of Genotyping




Wild type: ①PCR reaction obtains a single WT band; ②PCR reaction obtains a single WT band.
Heterozygote: ①PCR reaction obtains a WT band and a Targeted band; ②PCR reaction obtains a WT band and a Targeted band.
Homozygote: ①PCR reaction obtains a single Targeted band; ②PCR reaction obtains a single Targeted band.
Note: The sizes of WT and Targeted band are shown below.

2. Primer Information

PCR No.	Primer No.	Primer Name	Sequence	Band Size
①(5'arm)	F1	JS08549-Homer1-5wt-tF1	GAATTTGTTTCTCCATAAGACGCCTC	WT: 253bp Targeted:357bp
	R1	JS08549-Homer1-5wt-tR1	TCTCCACATTAGCACACGATGG	
②(3'arm)	F2	JS08549-Homer1-3wt-tF1	ATTTCCCAGTGCCAGTCAGTGGT	WT: 324bp Targeted:425bp
	R2	JS08549-Homer1-3wt-tR1	TTGTCTGACTGCAAATCAATAAATCTG	

3. Gel Image & Conclusion

--



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃（-0.5℃/cycle）	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	



Reviewer:Ting Sun

Date:2020-5-12