

## Genotyping Report

|           |           |             |              |                    |             |
|-----------|-----------|-------------|--------------|--------------------|-------------|
| Strain ID | T018921   | Strain Type | CKO(Cas9)    | Genetic Background | C57BL/6JGpt |
| Designer  | Zifan Lin | Gene Name   | <i>Kcna3</i> |                    |             |

### 1. Strategy of Genotyping



Wild type: ①PCR reaction obtains a single WT band; ②PCR reaction obtains a single WT band.

Heterozygote: ①PCR reaction obtains a WT band and a Targeted band; ②PCR reaction obtains a WT band and a Targeted band.

Homozygote: ①PCR reaction obtains a single Targeted band; ②PCR reaction obtains a single Targeted band.


Note: The sizes of WT and Targeted band are shown below.

### 2. Primer Information

| PCR No.  | Primer No. | Primer Name           | Sequence                  | Band Size                  |
|----------|------------|-----------------------|---------------------------|----------------------------|
| ①(5'arm) | F1         | JS09215-Kcna3-5wt-tF1 | ATGGCTTGAGGATGTGGACTTCTC  | WT:354bp<br>Targeted:458bp |
|          | R1         | JS09215-Kcna3-5wt-tR1 | GCTCTAGTTTGCTGGATACTCCTGC |                            |
| ②(3'arm) | F2         | JS09215-Kcna3-3wt-tF1 | TGAACCAAAAGGGAAATGCCC     | WT:268bp<br>Targeted:369bp |
|          | R2         | JS09215-Kcna3-3wt-tR1 | GTCCGAATTAGTAACCTCCGCTGG  |                            |

### 3. Gel Image & Conclusion

|  |
|--|
|  |
|--|



#### 4. PCR Condition

| PCR Reaction Component |                                                 |      |             |
|------------------------|-------------------------------------------------|------|-------------|
| Seg.                   | reaction component                              |      | Volume (μl) |
| 1                      | 2 × Taq Master Mix , Dye Plus, (Vazyme P112-03) |      | 12.5        |
| 2                      | ddH2O                                           |      | 9.5         |
| 3                      | Primer A(10pmol/μl)                             |      | 1           |
| 4                      | Primer B(10pmol/μl)                             |      | 1           |
| 5                      | Template(≈100ng/μl)                             |      | 1           |
| PCR program            |                                                 |      |             |
| Seg.                   | Temp.                                           | Time | Cycle       |
| 1                      | 95℃                                             | 5min | 20×         |
| 2                      | 98℃                                             | 30s  |             |
| 3                      | 65℃ (-0.5℃/cycle)                               | 30s  |             |
| 4                      | 72℃                                             | 45s  |             |
| 5                      | 98℃                                             | 30s  | 20×         |
| 6                      | 55℃                                             | 30s  |             |
| 7                      | 72℃                                             | 45s  |             |
| 8                      | 72℃                                             | 5min |             |
| 9                      | 10℃                                             | hold |             |

