


Genotyping Report

Strain ID	T004068	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	<i>Lta4h</i>		

1. Strategy of Genotyping




2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
PCR①	T004068-F1	GTAGGAGTTGGGCTCTGTGTTTTTC	WT:1633bp KO:483bp
	T004068-R1	ACACTAGGTCAACGACTTGACC	
PCR②	T004068-F2	GTAGGAGTTGGGCTCTGTGTTTTTC	WT:642bp KO:0bp
	T004068-R2	ACTAAAGGCAGATGTATGACGG	

3. Gel Image

cttggcaaaagcatctgtgCATCGCATTGTCTGAGTAGGTGATAACTTCGTATAATGTATGCTATACGAAGTTATATG
TTCTTGCCCAAGGTCAGTTGGGGGCCCGATCCGagggaaaaggggttc//ggcaggaattgaa---+93bp-1244bp+1bp---
Cctggtaggttgagggtcaagtcgtgacc



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ (-0.5℃/cycle)	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2021-7-4



江苏集萃药康生物科技股份有限公司
GemPharmatech Co.,Ltd