


Genotyping Report

Strain ID	T028084	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	<i>Nrg4</i>		

1. Strategy of Genotyping




2. Primer Information

PCR No.	Primer No.	Primer Name	Sequence	Band Size
PCR①备用 GC:28%	F1	T028084-F1A	TTTGCTAATGGAACATTTGGCGC	WT:47863bp KO:639bp
	R1	T028084-R1A	AGTTTATCAAGAGTGATTGTTGCTTCAG	
PCR②备用	F2	T028084-F2A	AGGCTCAGTAGGAGACCCTATTTC	WT:412bp KO:0bp
	R2	T028084-R2A	TTGAAACAGGGTTCATTAGAGCCC	

3. Gel Image

aataaaattagcataagaggtagtaattca-----47224bp-----cctgggctctaataaacctgtttcaaaaa



Note: P: Positive control; WT: Wildtype control; B: Blank control (ddH₂O); M: DNA Ladder

4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH2O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	20×
2	98℃	30s	
3	65℃ （-0.5℃/cycle）	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	