

Brwd3 Cas9-CKO Strategy

Designer: Jia Yu

. Jia Iu

Reviewer: Xiaojing Li

Design Date: 2020-10-12

Project Overview



Project Name Brwd3


Project type Cas9-CKO

Strain background C57BL/6JGpt

Conditional Knockout strategy



This model will use CRISPR/Cas9 technology to edit the *Brwd3* gene. The schematic diagram is as follows:



Technical routes



- ➤ The *Brwd3* gene has 5 transcripts. According to the structure of *Brwd3* gene, exon5-exon6 of *Brwd3*-205(ENSMUST00000150434.7) transcript is recommended as the knockout region. The region contains 250bp coding sequence. Knock out the region will result in disruption of protein function.
- ➤ In this project we use CRISPR/Cas9 technology to modify *Brwd3* gene. The brief process is as follows:sgRNA was transcribed in vitro, donor vector was constructed.Cas9, sgRNA and Donor were microinjected into the fertilized eggs of C57BL/6JGpt mice.Fertilized eggs were transplanted to obtain positive F0 mice which were confirmed by PCR and sequencing. A stable F1 generation mouse model was obtained by mating positive F0 generation mice with C57BL/6JGpt mice.
- > The flox mice was knocked out after mating with mice expressing Cre recombinase, resulting in the loss of function of the target gene in specific tissues and cell types.

Notice



- > According to the existing MGI data,male chimeras hemizygous for a gene trapped allele exhibit short tail buds, microcephaly and, in some cases, embryonic growth retardation.
- > The *Brwd3* gene is located on the ChrX. If the knockout mice are crossed with other mice strains to obtain double gene positive homozygous mouse offspring, please avoid the two genes on the same chromosome.
- > This strategy is designed based on genetic information in existing databases. Due to the complexity of biological processes, all risk of loxp insertion on gene transcription, RNA splicing and protein translation cannot be predicted at existing technological level.

Gene information (NCBI)



Brwd3 bromodomain and WD repeat domain containing 3 [Mus musculus (house mouse)]

Gene ID: 382236, updated on 13-Mar-2020

Summary



Official Symbol Brwd3 provided by MGI

Official Full Name bromodomain and WD repeat domain containing 3 provided by MGI

Primary source MGI:MGI:3029414

See related Ensembl:ENSMUSG00000063663

Gene type protein coding
RefSeq status VALIDATED
Organism Mus musculus

Lineage Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia;

Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus

Also known as Brodl, D030064D06Rik, Gm10452, Gm596

Expression Broad expression in placenta adult (RPKM 3.9), CNS E11.5 (RPKM 1.8) and 20 other tissuesSee more

Orthologs <u>human</u> <u>all</u>


Transcript information (Ensembl)



The gene has 5 transcripts, all transcripts are shown below:


Name	Transcript ID	bp	Protein	Biotype	CCDS	UniProt	Flags
Brwd3-205	ENSMUST00000150434.7	11619	<u>1799aa</u>	Protein coding	CCDS41100	A2AHJ4	TSL:5 GENCODE basic APPRIS P1
Brwd3-202	ENSMUST00000101283.3	1066	<u>181aa</u>	Protein coding	-	A2AHJ4	TSL:1 GENCODE basic
Brwd3-201	ENSMUST00000041866.11	5839	<u>84aa</u>	Nonsense mediated decay	-	G3UZT0	TSL:5
Brwd3-204	ENSMUST00000144521.7	5740	<u>60aa</u>	Nonsense mediated decay	-	G3UXF9	TSL:5
Brwd3-203	ENSMUST00000134434.1	1654	No protein	Retained intron	-	-	TSL:1

The strategy is based on the design of *Brwd3-205* transcript, the transcription is shown below:




Genomic location distribution






Protein domain





Mouse phenotype description(MGI)





Phenotypes affected by the gene are marked in blue.Data quoted from MGI database(http://www.informatics.jax.org/).

According to the existing MGI data,male chimeras hemizygous for a gene trapped allele exhibit short tail buds, microcephaly and, in some cases, embryonic growth retardation.



If you have any questions, you are welcome to inquire.

Tel: 025-5864 1534





